

[the 600 series]



600



610



620



medical/laboratory seating & equipment



600



610



620

Blood Drawing Chair Specifications

600

- seat height: 21"
- seat dimension: 18"Wx17"D
- overall height: 36"
- overall width: 28.5"
- armrest adjustment height: 27"-36"
- overall depth: 21"
- weight capacity: 300lbs
- weight: 30lbs

Options:

- cabinet (not shown, uses 620 frame)

610

- seat height: 22"
- seat dimension: 18.5"Wx18.5"D
- overall height: 37.5"
- overall width: 28.5"
- armrest adjustment height: 27"-36"
- overall depth: 29"
- weight capacity: 300lbs
- COM requirements 2 yd.
- weight: 49lbs

Options:

- cabinet (not shown, uses 620 frame)

620

- seat height: 22"
- seat dimension: 24.5"Wx18.5"D
- overall height: 37.5"
- overall width: 34"
- armrest adjustment height: 27"-36"
- overall depth: 29"
- weight capacity: 450lbs
- COM requirements 2 yd.
- weight: 55lbs



medical/laboratory seating & equipment

P.O. Box 5981 High Point, NC 27262
Toll-Free: 1-888-636-2660 Fax: (336) 476-4026