

[the 860 series]



862



863
[polished chrome base]



862
[alum. base w/toecaps]



861
[black composite base]



Lab Seating Specifications

860 Series Lab Stool:

- single lever hand adjustment
- seat height adjustment: 23" - 33"
- seat diameter: 15.75"
- base diameter: 24"
- high density fire resistant foam
- 3/4" thick composite substrate with steel reinforcement plate
- adjustable seat back
- seat back dimensions: 15.75" x 9.5"
- 18" diameter foot rest
- standard with casters
- soft wheel casters or glides optional