

[the 980 series]



983



983
[polished chrome base]



982
[alum. base w/toecaps]



981
[black composite base]

Physician Seating Specifications 980 Series Exam Stool:

- 360 degree foot ring adjustment
- standard height adj.: 17" to 22"
- seat diameter: 15 3/4"
- base diameter: 24"
- high density fire resistant foam
- 3/4" thick composite substrate with steel reinforcement plate

Series Options

Back Rest Option



Back Rest: (specify "98X-B")

- adjustable back
- seat back dimensions:
15 3/4" x 9 1/2"

Height Adjustment:

- 19" to 26" - specify "7" (ex. 98X-7)

Soft Wheel Casters or Glides

intensa

medical & laboratory seating

P.O. Box 5981 High Point, NC 27262
Toll-Free: 1-888-636-2660 Fax: (336) 476-4026