

## PHLEBOTOMY SEATING

# 600 SERIES PHLEBOTOMY CHAIR

Upholstered Bariatric Blood Draw/Phlebotomy Chair with Padded Arm Rest  
Intensa offers a bariatric blood draw/phlebotomy chair designed with both the patient and phlebotomist in mind. The sensible and affordable design offers a small footprint and modern style while keeping patient comfort and well-being the highest priority.

### 620 Features:

- Powder coat frame finish - Foghorn
- Fully upholstered padded arm rests
- Height adjustable arm rests

 CF Stinson Pre-Approved Fabrics

 620 Product Info




### 620 Options:

	item no. †
Single Flip Arm	620-xx
Dual Flip Arms	620-DA-xx
Wall Saver Wireform	620 WALL SAVER-xx

options available at additional cost

### 620 Specifications:

	inch	cm		inch	cm		yd	m
Overall Width	34.00	86.36	Seat Width	24.50	62.23	COM Requirement <sup>1</sup>	2.5	2.29
Overall Depth	29.00	73.66	Seat Depth	18.50	46.99		lb	kg
Overall Height	37.50	95.25	Seat Height	22.00	55.88		Weight Capacity	800
Arm Height -low	27.00	62.23	Back Width	25.50	64.77	EVENLY DISTRIBUTED		
Arm Height -high	36.00	66.04	Back Height	17.00	43.18	Shipping Weight	55	24.95
			Powder Coat Finish	foghorn				

<sup>‡</sup> Ships: Ground at no additional charge


<sup>†</sup> -xx in the product number is the placeholder for vinyl choice

<sup>1</sup> Submit COM sample for approval



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

The most current and complete product information can be found on our website.  
Information contained herein is subject to change without notice.

 Color variations may occur due to monitor or printer settings.

Intensa, Inc. | 1810 S. Elm Street | High Point, NC 27260 | USA  
+1 336.884.4096 | [intensa@grahamfield.com](mailto:intensa@grahamfield.com) | [intensa.net](http://intensa.net)

Manufactured in the USA

