

[503 series exam lighting]



503MAG



503



medical/laboratory seating & equipment



503MAG



503



medical/laboratory seating & equipment

P.O. Box 5981 High Point, NC 27262
Toll-Free: 1-888-636-2660 Fax: (336) 884-4007

503 Series Exam Light Specifications

503

- patented dual shade system
- 5-leg cast iron base with casters
- chrome shaft height range: 33"-63"
- polyvinyl coated steel goose neck (28")
- 360° lamp adjustment
- UL approved 12' heavy duty cord with hospital grade plug
- uses standard incandescent bulb up to 100 watts or a 60 watt halogen bulb (not included)

503MAG

- includes a second poly-vinyl coated goose neck (30") with a 5" diameter 2x magnifier

Available Options

- available colors: black, beige, gray
- stationary base in place of standard casters available (502)
- bulb options: CCH60, Par 20 LED