

[the 990 series]



993



993
[polished chrome base]



992
[alum. base w/toecaps]



991
[black composite base]

Physician Seating Specifications 990 Series Exam Stool:

- 360 degree hand ring adjustment
- standard height adj.: 17" to 22"
- seat diameter: 15.75"
- base diameter: 24"
- high density fire resistant foam
- 3/4" thick composite substrate with steel reinforcement plate

Series Options

18" Foot Rest Option



Foot Rest: (specify "99X-FR")

- polished chrome finish
- foot rest diameter: 18"

Height Adjustment:

- 19" to 26" - specify "7" (ex. 99X-7)
- 23" to 33" - specify "10" (ex. 99X-10)

Soft Wheel Casters or Glides

intensa

medical & laboratory seating

P.O. Box 5981 High Point, NC 27262
Toll-Free: 1-888-636-2660 Fax: (336) 884-4007

www.intensa.net