




Medical/Laboratory
Furniture & Equipment

PRODUCT SPECIFICATION SHEET

860 Series Laboratory Seating

861



All 860 series laboratory stools offer Pneumatic Height Adjustment, Lateral Back Rest Adjustment and can safely support 350 lbs.



862



863

860 Series Specifications

Lab Seating Chair

3/4" Thick Composite Substrate with Steel Reinforcement Plate

Seat Diameter: 16"

Seat Height: 23" - 33"

Seat Foam Thickness: 3"

Back Dimensions: 15.75" H x 9.75" W

Base Diameter: 24"

20" Diameter Foot Ring

860 Series Standard Features & Available Options

Adjustable Seat Back

Single Lever Height Adjustment

High-Density Fire-Resistant Foam

Standard with Hard Wheel Casters

COM Requirements: 1 yd.

Available Options:

California TB-133

Soft Wheel Casters or Glides

Crescent Back (-CB)